

DRAGER NARKOMED 2B




Certified Refurbished

The smart choice for reliable anesthesia

The NAD Narkomed 2B is one of Drager's longest running and most reliable anesthesia machines. It comes standard with O₂ & N₂O and holds three vaporizers. The Narkomed 2B also includes F-Style ascending bellows with oxygen, volume and pressure monitoring.

► A Quick Look

- 19.1 and 19.3 vaporizers available
- ORM standard; ORMC optional
- Ventilator options: Pneumatic or electronic, ascending or descending bellows
- Oxygen monitor and CO₂ absorber standard
- AVE Electronic ventilator
- Absorber and scavenger
- Oxygen analyzer
- Breathing pressure monitor
- Respiratory volume
- Centralized alarm and data
- Oxygen ratio controller
- Pressure limit controller



Drawer pack configurations may vary from picture.

Technical Specifications

Maximum Height	68" (173 cm)
Maximum Width	40" (102 cm)
Maximum Depth	25" (63.5 cm)
Table-top Height	34" (86 cm)
Table-top Area	271 in. sq. (1,749 cm ²)
Shelf Area	200 in. sq. (1,290 cm ²)
Drawer Dimensions	
– Large	14"x21"x7" (36 x 52 x 18 cm)
– Small	14"x17"x3" (36 x 43 x 8 cm)
Nominal Weight	325 lb. (147 kg)



IPH MAG
www.magbiomedica.com