

OHMEDA EXCEL 110

Product Description

An Excel 110 is a smaller version of the 210 models, without medical air. Ohmeda 7000 or 7800 ventilators and Tec series vaporizers are optional on the Excel 110, while the GMS CO2 Absorber and Hypoxic Guard are standard. Drawer pack configurations may vary from picture.

Features:

- Compact Frame; Fits Space-Constrained Areas.
- Integrated Oxygen, Volume, Airway Pressure Monitoring.
- Enclosed, Non-interchangeable Pneumatic Circuitry.
- Differentiated Pipeline, Cylinder Gauges.
- Master On / Off Switch; Activates Electronics and Pneumatics.
- Link-25 Hypoxic Mixture Guard System.
- Pin-Indexed, Gas-Specific Dual O2 and N2O Flowmeters.
- Selectatec® Vaporizer Manifold for two Tech 5 Vaporizers.
- Guarded O2 Flush Button.
- Oxygen Power Outlet.
- Fixed shelving with added Monitoring space.
- Two-drawer Storage Cabinet.
- Single-action Brake / Footrest.
- Full-height Dovetail Accessory Mounting System.
- Large, Easy-running Casters.

Specifications:

- H x W x D - 137 x 66 x 66 cm / 54 x 26 x 26 in
- Weight – 122 kg / 270 lbs
- Height – 168 cm / 66 in
- Width of Frame – 50.7 cm / 20 in
- Middle Shelf Maximum Load - 11 kg / 24 lbs
- Middle Shelf - 42 x 31.7 cm / 16.6 x 12.5 in



IPHMAG
www.magbiomedica.com