

DRAGER NARKOMED 4

The North American Drager Narkomed 4 anesthesia system is a very popular system for anesthesia delivery, monitoring and management. At the time of this posting the Drager Narkomed 4 is one of the most advanced NAD machine available, and is equipped for enhanced critical event avoidance and rapid response.

Standard features include all NK2B features plus Touch Screen, Remote Display, SpO2, non-invasive blood pressure, EtCO2 and agent monitoring.

Features:

- 19.1 vaporizers
- AVE ventilator
- NIB
- SPO2
- CO2
- Built in monitor.

Options:

- Air system
- PEEP
- Boom Arm.

Specifications:

- Weight: 500lb (226.8kg)
- Height: 64 3/4 in (165.5 cm)
- Width: 31 3/4 in (80.6 cm)
- Depth: 25 1/8 in (63.8 cm)
- Shelf area: 200 in sq (1,290 cm sq)
- Power Supply: 117 volts
- Optional Power Supply: 220/240 volts
- Backup Battery:
 - Charging time: <12 hours
 - Reserve power time: >25 minutes



IPHMAG
www.magbiomedica.com