

Bed 850-852

Refurbished Hill-Rom Centra 850/852 hospital beds are adjustable and have a sturdy, long-lasting frame.

Perfect for hospitals, surgery centers and teaching facilities!

Specifications

- Three motors
- High-low travel time: 21 seconds
- Sleep surface: 36" x 82"
- Bed frame clearance with 5" casters: 7.5"
- Low position: 17"
- High position: 30"
- Central pedal for braking and steering
- Night light
- CPR Release

Options

- Radio/TV
- Bed exit
- Nurse call
- Parallelogram (new style) side rails
- Orthopedic (hard pan) surface
- Mattress (built to OEM specifications)
- Light neutral paint
- Roller bumpers

Refurbishing Procedures

- Spring sleep surface (standard) or orthopedic/hard pan surface (upgrade)
- New PC boards with a 5-year warranty
- Fully rebuilt side rails w/ new plastics and switches
- New head/foot boards, your choice of laminate colors, lifetime warranty
- New paint: Taupe (standard) or light neutral (upgrade)
- New labels
- New motor couplings and capacitors
- New mattresses built to OEM specifications
- Rebuilt casters with new components and wheels
- 220-Volt transformers available
- Full function tests
- Basic bed function (head, knee, high-low)



IPHMAG

www.magbiomedica.com