

# Air shields isolette C550 QT-XL



Certified Refurbished

## Easy-to-use infant incubator with a simple, easy-to-read display

- Provides stable temperature
- Mattress tray tilt mechanism
- Variable height adjustment
- Baby-temperature and air-temperature controls



### Specifications

**Dimensions** ..... 140.4 mm x 119.4 mm x 59.4 mm

**Weight** ..... 140.6 kg

**Power Supply** ..... 240 V 50 Hz ac

### Access

- Hand ports ..... 5
- Doors ..... 1 (front)
- Tubing ports ..... 6
- Tray withdrawal ..... Front
- Tray tilting mechanism ..... Pneumatic  $\pm 9^\circ$

### Controller

- Air-temperature ..... 20° C to 39° C (0.1° C)
- Baby-temperature ..... 34° C to 38° C (0.1° C)

### Alarms

- Mains power failure
- High air temperature
- Skin temperature

IPHMAG

[www.magbiomedica.com](http://www.magbiomedica.com)