

CODEMASTER XL



Certified Refurbished

Delivers the high-speed and high-performance needed in general care

► A Quick Look

- High performance defibrillator
- Cardioscope and recorder
- Available with a noninvasive pacemaker
- Rapid charge
- Light weight
- True Easy 1-2-3 Operation
- Built-in pediatric paddles
- Adjustable heart rate alarms
- High capacity battery
- Patient contact indicator
- Includes all patient cables, a new battery and operator's manual

The Codemaster XL+ combines a DC defibrillator/ cardioscope, an annotated recorder and hands-off defibrillation capabilities into a single, high-performance unit. It operates on AC power or internal battery supply. It is also available with a noninvasive pacemaker.

Designed for high-use areas of the hospital, the high-speed Codemaster XL+ is perfect for operating rooms, intensive care units, cardiac care units & emergency departments. It charges from zero to 360 joules in under five seconds and offers a variety of advanced options so you can easily customize it to meet your facility's needs.

Order the patient-ready Codemaster XL+ from DRE and receive a one year defective parts warranty.

IPHMAG

www.magbiomedica.com