

VALLEYLAB FORCE2



Certified Refurbished



2 surgeons
can work
simultaneously
from a single generator
with the Force 2

The Force 2 ESU electrocautery unit provides cutting and coagulation in monopolar or bipolar outputs.

- Audible and visual alarms, including REM™ Fault Indicator.
- REM™ (Return Electrode Monitoring) comes standard.
- Available with 300 watts for output in the cut setting.
- Includes a monopolar footswitch (a bipolar foot switch is optional)

► A Quick Look

Blended cut modes

Provide flexibility through varying degrees of hemostasis. Choose one of three preset modes.

Simultaneous independent coagulation

Permits two surgeons to fulgurate from a single generator, for added convenience and efficiency. Accessories are activated only when keyed, reducing the potential for injury caused by inadvertent activation.

Low voltage coagulation

Ensures controlled, precise desiccation with less destruction of peripheral tissue, making the Force 2 generator ideal for laparoscopic procedures.

WEIGHT AND DIMENSIONS

Height: 8" (20 cm) / Width: 13" (33 cm)
Length: 21" (53 cm) / Weight: 23 lbs. (10.4 kg)



Available with an optional ESU cart

IPHMAG
www.magbiomedica.com