



For light
that's as

BRIGHT

as you need

► **A Quick Look**

- Four settings of electronic dimming with outputs of 100%, 85%, 70% and 50%
- Automatically switches to reserve bulb in case of primary bulb failure
- Rotational contacts at key pivots
- Available in single, dual or triple ceiling mounted

The DRE Maxx Luxx II provides four brightness settings, allowing you to adjust its output of cool, color-corrected light to fit your needs. The DRE Maxx Luxx II also includes adjustable pattern sizes; electronic touch controls for on/off and dimming; smooth, flexible and extensive suspension arm for maximum coverage; a reserve bulb that automatically switches over in the event of a primary bulb failure; and excellent shadow control.

The Maxx Luxx II Video Camera

features a 27x optical zoom and 10x digital zoom with auto focus, manual focus, and the option of a single-input controller or one that accommodates up to 8 video cameras.

3^{ceiling mounting} options available

single, double and triple



Product Specifications

Foot Candles	9,500 / light head at 1m	Rated Bulb Life (primary and secondary)	1000 hours
Lux	102,000 @ 1 Meter	Light of Dimming Settings (discrete)	100%, 85%, 70%, 50%
Kelvin	4000° color temp	Reflector	20" Facetted Alzak Alum.
Depth of Field	48" (1219mm)	Filter Type	I/R Filter Glass
Focus Distance	1,400 mm	Sterilizable Handle	Standard
Light Pattern	165mm to 229mm	Input Voltage Range	100-125 Vac 50/60 Hz
Extended arm coverage from center			220-240 Vac 50/60 Hz
of mount to center of light head	74" (1880 mm)	Lamp Voltage	23-25 Vac
Rotation at ceiling pivot	540 degrees	Light Head Dimensions	23" Dia. x 14" High
Minimum mounting height	96" (2438 mm)	Arm Length (max length)	59"
Maximum mounting height	156" (3962 mm)	Lamp Head Rotation	Continuous
Multi-faceted reflector	20" diameter (508 mm)	Vertical Arm Movement	+20, -70 degrees
Outer head	23" diameter (584 mm)	Ceiling Mount/Extension Arm	Continuous
Light Intensity at 39"	102,000 lux (9500 ft. candles)	Articulating/Extension Arm Joint	Continuous
Light Field Diameter (adjustable)	6.5"-9.0"	Articulating Arm/Yoke Joint	Continuous
Depth of Illumination	48"	Yoke/Light Head Joint	± 150°
Color of Temperature	4000 K	Weights (based on 8 ft. ceiling height):	
Color of Rendering Index	90 CRI	Single Ceiling Mount (w/ arms)	36 lbs.
Light Source (primary)	150 W Tungsten Halogen Bulb	Double Ceiling Mount (w/ arms)	62 lbs.
Light Source		Trio Ceiling Mount (w/ arms)	88 lbs.
(secondary, back-up)	150 W Tungsten Halogen Bulb	Light Head Assembly	17.5 lbs.