


**Quality, Reliability, and
Innovation**

Surgical Table Systems
from STERIS Corporation



STERIS®




Cmax™ Surgical Table

Precise mobility and positioning capabilities

- > Radiolucent, fully posturing tabletop
- > 18" (45 cm) longitudinal slide
- > Intelligent hand control with digital communication display
- > 1100 lb (500 kg) patient support, including raise/lower (table top not articulated or slid)
- > 1000 lb (454 kg) full patient posturing (table top not slid)
- > 600 lb (272 kg) full patient posturing, including table top slide



Amsco® 3085 SP Surgical Table

The ultimate surgical platform

- > 1000 lb. (454 kg) patient weight capacity
- > Radiolucent, fully posturing tabletop
- > The industry standard in versatility & reliability
- > Accessories for every procedure


STERIS surgical tables, accessories, and support systems provide complete flexibility for patient positioning, no matter how difficult the procedural or imaging requirements. These products assure you increased productivity, ease of use, and long term reliability. No other surgical table system better meets the challenges of today's operating room.



Amsco® OrthoVision® Surgical Table

The industry standard for orthopedic tables

- > Excellent orthopedic positioning capabilities
- > Radiolucent abductor bars allow for quality imaging
- > Exclusive translating top
- > 400 lb. (181 kg) patient weight capacity



SurgiGraphic® 6000 Image Guided Surgical Table

A breakthrough for image guided surgery

- > Complete surgeon control for precise positioning
- > Radiology-quality imaging for the O.R.
- > For endovascular, pain management, and other image guided procedures
- > 500 lb. (227 kg) patient weight capacity



SurgiGraphic® 1027 Image Guided Surgical Table

Total C-arm access, uncompromising value

- > Portable imaging capability
- > Open design allows 360° C-arm access
- > For endovascular, pain management, and other image guided procedures
- > 500 lb. (227 kg) patient weight capacity

STERIS: A complete line of surgical table solutions



Cmax™ Surgical Table



Amsco® 3085 SP Surgical Table and Orthopedic Extension



Amsco® OrthoVision® Surgical Table



SurgiGraphic® 6000 Image Guided Surgical Table



SurgiGraphic® 1027 Image Guided Surgical Table

Specifications

	Cmax table	Amsco 3085 SP table	Amsco OrthoVision table	SurgiGraphic 6000 table	SurgiGraphic 1027 table
Table top length	81" (206 cm)	76" (193 cm)	31" (79 cm)	92" (234 cm)	90" (229 cm)
Table top width	22" (55 cm)	20" (51 cm)	20" tapered (51 cm)	21" (53 cm)	21" (53 cm)
Table height range Minimum/Maximum	26" (68 cm) 45" (114 cm)	27" (69 cm) / 44" (112 cm)	29" (74 cm) / 47" (120 cm)	32" (81 cm) / 42" (107 cm)	27" (69 cm) / 40" (102 cm)
Trendelenburg / Reverse Trendelenburg	30°	25°	15°	16°	12°
Lateral tilt (left & right)	20°	18°	15°	16°	18°
Patient weight capacity	1100 lb. (500 kg) (incl. raise/lower)	1000 lb. (454 kg)	400 lb. (181 kg)	500 lb. (227 kg)	500 lb. (227 kg)
Sliding tabletop	18" (45 cm) longitudinal side	N/A	4" (10 cm) lateral translation	16" (41 cm) longitudinal float 8" (20 cm) lateral float	N/A

STERIS OFFICES WORLDWIDE

Benelux	32 2 523 2488	Japan	81 78 321 2271
Canada	800 661 3937	Korea	82 2 517 1517
China	86 21 6137 1166	Latin America	800 884 9550
France	33 2 38 70 83 50	Malaysia	6 03 7954 9822
Germany	49 2203 890 6969	Singapore	65 68 41 7677
Greece	30 210 6800848	Spain	34 916 585 920
India	91 33 2367 5150	Switzerland	32 376 02 0200
Italy	39 022 130341	United Kingdom	44 1256 840400

Document # #M2695EN.2009-01, Rev. C
Printed 02/2009, 500

©2009 STERIS Corporation.
All rights reserved. Printed in USA.

STERIS®



STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 ■ USA
440-354-2600 ■ 800-548-4873
www.steris.com