

INTRODUCING

STREAMLINE

X-Ray imaging tables from Surgical Tables Inc.

DESIGN

This table line is designed for high patient throughput for surgery centers, pain management, and other multidisciplinary medical facilities.


VALUE

The Streamline series is the lowest priced table in the imaging table industry. This line along with the best selling Max series and EconoMax series tables, makes STI the company with the most table options to choose from in the medical table business.



Carbon fiber diving board style imaging table

"absolutely the best price and value in the table industry"



FEATURES

- Diving board style for easy c-arm access
- Large patient capacity – 500 pound load
- 2" patient high quality comfort pad
- Facial Cut-out for prone procedures
- Light weight - easy to physically move the table
- 24" top eliminates need for arm boards
- Simple to use top quality wheel locks
- Easy patient transfer to table




Tel: 888-737-5044

Tel: 978-777-4031

Fax: 978-750-0877

3 Ajoon Way Unit D10
Middleton, MA 01949

www.surgicaltablesinc.com



Carbon fiber top of 24" width eliminates need for armboards and 84" length can accommodate even the tallest patients

Facial Access Hole

2" thick high quality patient comfort pad

Easy to use handswitch with power on indicator

High end actuator motor can lift over 500 pound load

Diving board style easy c-arm access



Lightweight which makes it easy for staff to move the table around the room or room to room



For a price quote or for more info call

888-737-5044

Or check out our website at www.surgicaltablesinc.com

OPTIONS ▾

MODULAR SIDE RAIL SYSTEM (Removeable armboard rotates 180 degrees)



SURGICAL ARMBORDS (Removeable armboard rotates 180 degrees)



PROCEDURAL ARMBORDS (Each armboard rotates 180 degrees)



CERVICAL IMAGING PLATFORM



FOOT CONTROLS

